



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat MMS19 Monoclonal Antibody

单克隆抗体

产品货号: MA02039

第三版

免疫原:Fusion protein of MMS19

亚型:IgG2b

克隆:12

**描述:**MMS19 (MMS19 nucleotide excision repair homolog),also known as MET18,is a 1,030 amino acid nuclear protein containing seven HEAT repeats that belongs to the MET18/MMS19 family. Via its interactions with TFIIH p80 and TFIIH p89 helicases,MMS19 plays a role in nucleotide excision repair (NER) and RNA polymerase II (Pol II) transcription. MMS19 may also function as a transcriptional coactivator of estrogen receptor. While ubiquitously expressed,highest levels of MMS19 have been found in testis. At least five distinct MMS19 protein isoforms exist,which are produced by alternative splicing events.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF