

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MSH6 Monoclonal Antibody

单克隆抗体

产品货号: MA02042

第三版

免疫原:Fusion protein of MSH6

亚型:IgG2a

克隆:120

描述:MSH6,also named as DNA mismatch repair protein Msh6 or GTMBP,is a 1360 amino acid protein,which contains one PWWP domain and belongs to the DNA mismatch repair MutS family. MSH6 localizes in the nucleus and is a component of the post-replicative DNA mismatch repair system. Msh2 and Msh6 form a protein complex required to repair mismatches generated during DNA replication.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF