

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse MTA2 Monoclonal Antibody

单克隆抗体

产品货号: MA02043

第三版

免疫原:Fusion protein of MTA2

亚型:IgG2b

克隆:6101

描述:Metastasis-associated protein MTA2,also named as MTA1L1 or PID,is a 668 amino acid protein,which contains one BAH domain,one ELM2 domain,one GATA-type zinc finger and one SANT domain. MTA2 localizes in the nucleus and is widely expressed in many tissues. MTA2 may be involved in the regulation of gene expression as repressor and activator.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF