



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human MYH2 Monoclonal Antibody

单克隆抗体

产品货号: MA02046

第三版

免疫原:Synthetic Peptide

亚型:IgG2b

克隆:47

描述:MYH2 (Myosin II) is a member of the class II or conventional myosin heavy chains. Myosin II was first isolated from muscle but is also found in non muscle cells, and it is especially enriched in highly motile cell types such as amoebae. It forms bipolar filaments that interact with actin filaments to produce contraction. And it is the motor protein that generates force to drive muscle contraction. Functions in skeletal muscle contraction. The antibody is specific to MYH2, will not bind other myosins.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF