



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NAPRT Monoclonal Antibody

单克隆抗体

产品货号: MA02048

第三版

免疫原:Fusion protein of NAPRT1

亚型:IgG2a

克隆:77

描述:Nicotinic acid (NA) is a coenzyme in cellular redox reactions,and is an essential component of metabolic pathways in all living cells. NAPRT1 (Nicotinate phosphoribosyltransferase) is essential for increasing cellular NAD levels and,thus,to prevent oxidative stress of cells. NAPRT1 converts Nicotinic acid (NA; niacin) to NA mononucleotide (NaMN),which is then converted to NA adenine dinucleotide (NaAD),and finally to nicotinamide adenine dinucleotide (NAD).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF