



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse NAPSA Monoclonal Antibody

单克隆抗体

产品货号: MA02049

第三版

免疫原:Fusion protein of Napsin

亚型:IgG1

克隆:779

**描述:**Napsin is found in two isoforms,napsin A and B,with highly homologous nucleotide sequences (91.2%). Napsin A appears to be a functional proteinase,predominantly expressed in lung and kidney. Napsin B is transcribed exclusively in cells related to the immune system and lacks an in-frame stop codon and is believed to be a pseudogene.. Napsin A is superior to TTF-1 in distinguishing primary lung ACA from other carcinomas (except kidney),particularly primary lung small cell carcinoma,and primary thyroid carcinoma..

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF