

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NCALD Monoclonal Antibody

单克隆抗体

产品货号: MA02051

第三版

免疫原:Fusion protein of NCALD

亚型:IgG2a

克隆:213

描述:NCALD(Neurocalcin-delta) is a member of the EF-hand calcium-binding protein superfamily.It may be involved in the calcium-dependent regulation of rhodopsin phosphorylation.It is a high-affinityCa²⁺ sensor (Kd of approx. 0.6 uM),which is most abundantly expressed in the Purkinje cells of the cerebellum.This antibody is specific to NCALD.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF