



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat LXR alpha Monoclonal Antibody

单克隆抗体

产品货号: MA02059

第三版

免疫原:Fusion protein of NR1H3

亚型:IgG1

克隆:31

描述:NR1H3,also known as LXRA,is a liver X receptors,which form a subfamily of the nuclear receptor superfamily and are key regulators of macrophage function,controlling transcriptional programs involved in lipid homeostasis and inflammation.. It's also a ligand-activated LXRs and forms heterodimers with retinoid X receptor,then participates in regulating expression of target genes containing LXR response elements.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF