
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ENO2 Monoclonal Antibody

单克隆抗体

产品货号: MA02062

第三版

免疫原:Fusion protein of NSE

亚型:IgG1

克隆:345

描述:Enolases are cytoplasmic glycolytic enzymes that may be involved in differentiation. Three enolase isoenzymes, alpha, beta and gamma, are present in mammalian cells. The alpha form is expressed in most tissues, whereas the beta form is expressed in muscle tissue. The gamma enolase (ENO2), a homodimer, is found in mature neurons and cells of neuronal origin. A switch from alpha enolase to gamma enolase occurs in neural tissue during development.

配方:

浓度:None

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF