



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human NUMBL Monoclonal Antibody

单克隆抗体

产品货号: MA02064

第三版

免疫原:Fusion protein of NUMBL

亚型:IgG2a

克隆:301

描述:Numbl (also known as numbl like or Nbl), as a conserved homolog of Drosophila Numb, is specifically expressed in the brain. Numbl is a cytoplasm protein and has redundant functions in embryonic neurogenesis. It plays an important role in ependymal wall integrity and subventricular zone neuroblast survival. However, its expression and function in the central nervous system lesion are still unclear.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF