

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ORA11 Monoclonal Antibody

单克隆抗体

产品货号: MA02067

第三版

免疫原:Fusion protein of ORA11

亚型:IgG2a

克隆:8998

描述:Store-operated calcium channels (SOCs) play essential functions in various physiological processes. Orai1 (also named CRACM1) is a critical molecular component of SOCs. Orai1 protein is a 32-55 kDa plasma membrane protein with 4 transmembrane domains. Human patients with a severe combined immunodeficiency (SCID) syndrome possess point mutations in Orai1 gene.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF