

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ORM1 Monoclonal Antibody

单克隆抗体

产品货号: MA02068

第三版

免疫原:Fusion protein of ORM1/2

亚型:IgG1

克隆:321

**描述:**Alpha-1-acid glycoprotein 1 (AGP1),also called orosomucoid-1 (ORM1),is a glycoprotein synthesized mostly by hepatocytes and present in human plasma. ORM1 is an acute-phase reactant protein controlled by glucocorticoids,interleukin-1 and interleukin-6,and increase up to 5-50 times upon infection and/or inflammation. Anti-apoptotic effect and role as immunomodulator of ORM have been reported. ORM is an important carrier for synthetic drugs and influences their distribution and availability in the body. This antibody recognizes a band about 44 kDa in human plasma which may be due to the glycosylation of ORM1 or the dimer formation of the protein. This antibody recognizes both ORM1 and ORM2.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF