

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human ORM2 Monoclonal Antibody

单克隆抗体

产品货号: MA02069

第三版

免疫原:Fusion protein of ORM2-Specific

亚型:IgG1

克隆:495

描述:ORM2,also named as AGP2 and OMD 2,belongs to the calycin superfamily and Lipocalin family. It functions as transport protein in the blood stream. ORM2 binds various hydrophobic ligands in the interior of its beta-barrel domain. It also binds synthetic drugs and influences their distribution and availability. ORM2 appears to function in modulating the activity of the immune system during the acute-phase reaction. ORM2 is one of the conserved endoplasmic reticulum membrane proteins which regulating lipid homeostasis and protein quality control. This antibody can specially recognize ORM2. It recognizes a band about 45 kDa in human plasma which may be due to the glycosylation of ORM2 or the dimer formation of the protein.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF