

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse/Rat P53 Monoclonal Antibody

单克隆抗体

产品货号: MA02074

---

第三版

免疫原:Fusion protein of P53

亚型:IgG2b

克隆:665

描述:TP53,also named as P53 and NY-CO-13,belongs to the p53 family. It has 9 isoforms. In SDS-Page,the MW is about 53kd. TP53 acts as a tumor suppressor in many tumor types; induces growth arrest or apoptosis depending on the physiological circumstances and cell type. It is involved in cell cycle regulation as a trans-activator that acts to negatively regulate cell division by controlling a set of genes required for this process. The antibody reacts specifically with 53 kDa protein.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF