



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat Pan Acetylation Monoclonal Antibody

单克隆抗体

产品货号: MA02075

第三版

免疫原:Synthetic Peptide

亚型:IgG1

克隆:774

描述:Acetylation and De-acetylation are the important reaction in epigenetic. Many proteins,especially Histone,are with Acetylated of Lysine. This antibody is used to detect the Acetylation of Lysine. The ~17 kDa in our detection is Acetylated Histone.

配方:

浓度:1 mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF