



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse PAX1 Monoclonal Antibody

单克隆抗体

产品货号: MA02076

第三版

免疫原:Synthetic Peptide

亚型:IgG2b

克隆:302

描述:PAX1,also named as HUP48,is a transcriptional activator. It may play a role in the formation of segmented structures of the embryo and may play an important role in the normal development of the vertebral column. It is a marker of the sclerotomal compartment. PAX1 was reported that,quantitative measurement of PAX1 hyper-methylation in cervical scrapings is highly sensitive and is more specific than HC2 in detection of cervical cancer. It is expressed by vSMC and pericytes during the earliest phases of vascular wall formation. This antibody is specific to PAX1.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF