



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PAX8 Monoclonal Antibody

单克隆抗体

产品货号: MA02078

第三版

免疫原:Synthetic Peptide

亚型:IgG2b

克隆:344

描述:PAX8 is a member of the paired box (PAX) family of transcription factors, typically containing a paired box domain and a paired-type homeodomain. The specific function of the PAX8 is unknown but it may involve kidney cell differentiation, thyroid development, or thyroid dysgenesis. Several other splice variants have been proposed but the full nature of these products has not been described. Pax8 is also a marker of otic progenitor cells. Mutation of PAX8 will cause congenital hypothyroidism non-goitrous type 2 (CHNG2). This antibody is specific to Pax8. It will not cross react with other PAX members

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF