



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PDCD4 Monoclonal Antibody

单克隆抗体

产品货号: MA02080

第三版

免疫原:Fusion protein of PDCD4

亚型:IgG1

克隆:911

描述:Programmed cell death 4 (Pdc4) is a novel tumor suppressor that inhibits translation, progression and invasion. It was first identified as being differentially upregulated during apoptosis. Pdc4 interferes with the activity of the eukaryotic initiation factor (eIF) 4A by displacing the scaffold protein eIF4G from its binding to the RNA helicase eIF4A.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF