



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PDHA1 Monoclonal Antibody

单克隆抗体

产品货号: MA02081

第三版

免疫原:Fusion protein of PDHA1

亚型:IgG2a

克隆:435

描述:PDHA1(Pyruvate dehydrogenase E1 component subunit alpha,somatic form,mitochondrial) is also named as PHE1A.It is one of the 3 enzymes of the pyruvate dehydrogenase complex which is a nuclear-encoded mitochondrial matrix multienzyme complex that provides the primary link between glycolysis and the tricarboxylic acid (TCA) cycle by catalyzing the irreversible conversion of pyruvate into acetyl-CoA.It has 4 isoforms produced by alternative splicing.Defects in PDHA1 are a cause of pyruvate dehydrogenase E1-alpha deficiency (PDHAD) and X-linked Leigh syndrome (X-LS).

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF