

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PFN2 Monoclonal Antibody

单克隆抗体

产品货号: MA02084

第三版

免疫原:Fusion protein of PFN2

亚型:IgG2a

克隆:201

描述:PFN2 (profilin 2) is a ubiquitous actin monomer-binding protein belonging to the profilin family. It is thought to regulate actin polymerization in response to extracellular signals. The MW of this protein is 15 kDa, and this antibody specially recognises the 15 kDa protein.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF