

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PGM1 Monoclonal Antibody

单克隆抗体

产品货号: MA02085

第三版

免疫原:Fusion protein of PGM1

亚型:IgG2b

克隆:312

描述:PGM1(Phosphoglucomutase-1) is also named as glucose phosphomutase 1 and belongs to the phosphohexose mutase family. It catalyzes the transfer of phosphate between the 1 and 6 positions of glucose. In most cell types,PGM1 isozymes predominate,representing about 90% of total PGM activity. One exception is red cells,where PGM2 is a major isozyme. Defects in PGM1 are the cause of glycogen storage disease type 14 (GSD14). It has 2 isoforms produced by alternative splicing.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF