

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat PLAA Monoclonal Antibody

单克隆抗体

产品货号: MA02090

第三版

免疫原:Fusion protein of PLAA

亚型:IgG1

克隆:555

**描述:**Phospholipase A-2-activating protein (PLAA),also named PLAP,can interact with ubiquitin and is involved in the maintenance of ubiquitin levels. PLAA contains multiple consensus domains including WD repeats,PFU domain and PUL domain. Although the function of PLAA remains unclear,ubiquitin binding of PLAA might be the central role that connects ubiquitination and degradation in ERAD. Recent finding revealed that PLAA served as a novel nociveptive mediator after incision,and its expression level is regulated by miR-203.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF