



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human Placental Alkaline Phosphatase Monoclonal Antibody

单克隆抗体

产品货号: MA02091

第三版

免疫原:Fusion protein of PLAP

亚型:IgG1

克隆:564

描述:Placental alkaline phosphatase (PLAP) is a glycosyl phosphatidylinositol (GPI)-anchored sialoglycoprotein ,which shows greater resistance to heat inactivation than other alkaline phosphatase. Presented at high levels in placental trophoblasts,PLAP is recognized as an important marker of differentiation in human malignancies.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF