

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PLOD3 Monoclonal Antibody

单克隆抗体

产品货号: MA02092

第三版

免疫原:Fusion protein of PLOD3

亚型:IgG2a

克隆:367

描述:PLOD3,also named as LH3,forms hydroxylysine residues in -Xaa-Lys-Gly- sequences in collagens. These hydroxylysines serve as sites of attachment for carbohydrate units and are essential for the stability of the intermolecular collagen cross-links. The major function of PLOD3 in osteoblasts is to glucosylate galactosylhydroxylysine residues in type I collagen.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF