



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PGR Monoclonal Antibody

单克隆抗体

产品货号: MA02094

第三版

免疫原:Fusion protein of PR

亚型:IgG1

克隆:227

描述:PGR,also named as NR3C3 and PR,belongs to the nuclear hormone receptor family and NR3 subfamily. The steroid hormones and their receptors are involved in the regulation of eukaryotic gene expression and affect cellular proliferation and differentiation in target tissues. Progesterone receptor isoform B (PRB) is involved activation of c-SRC/MAPK signaling on hormone stimulation. The ER and PR status has been used as a predictor of breast carcinoma responsiveness to endocrine therapy and as a prognostic indicator for early recurrence.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF