

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PRDX4 Monoclonal Antibody

单克隆抗体

产品货号: MA02095

第三版

免疫原:Fusion protein of PRDX4

亚型:IgG2b

克隆:473

描述:PRDX4(Peroxiredoxin-4) is also named as AOE37-2,Prx-IV and belongs to the AhpC/TSA family. PRDX4 is associated with acrosome formation during rat spermatogenesis and has a protective role in the male reproductive tract,because the phenotype of mice lacking this isoform includes testicular atrophy and increased sperm DNA damage. It is initially synthesized as a membrane-binding 31-kDa protein and processed into a 27-kDa secretory form and is discarded with the residual bodies. PRDX4 is a pentamer of dimers.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF