



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human IL18 Monoclonal Antibody

单克隆抗体

产品货号: MA02097

第三版

免疫原:Fusion protein of Pre-IL18

亚型:IgG2b

克隆:112

描述:IL18, is a proinflammatory cytokine involved in the development of Th1 cells and in immune response. It can stimulate the NK cells and certain T cells to release interferon gamma which plays an important role in activating the macrophages or other cells. IL18 has been demonstrated to have the potential to enhance Fas ligand-mediated cytotoxicity, which is increased in PE and regulates placental apoptosis. IL18 is synthesized as a 24 kDa precursor and then cleaved into a biologically active 18 kDa form.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF