



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Rat PHB Monoclonal Antibody

单克隆抗体

产品货号: MA02098

第三版

免疫原:Fusion protein of Prohibitin

亚型:IgG1

克隆:73

描述:Prohibitin,also known as PHB,a member of the Band-7 family of proteins,is widely distributed in bacteria,plants,yeast,protozoa and mammals and is important in cell proliferation,differentiation and apoptosis . This protein localizes to the inner membrane of mitochondria,where it acts as a chaperone protein,and is also found in the nucleus,where it negatively regulates transcription . Recent study confirmed that the expression and distribution of PHB,which is a nuclear matrix protein,affect the apoptosis of HaCaT cells and its co-localization with specific gene products connected with cell apoptosis .

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF