



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SERPINA5 Monoclonal Antibody

单克隆抗体

产品货号: MA02099

第三版

免疫原:Fusion protein of Protein C inhibitor

亚型:IgG2a

克隆:34

描述:SERPINA5 (also called Protein C Inhibitor,PCI ) is a serpin type of serine protease inhibitor that is found in many tissues and fluids in human,including blood plasma,seminal plasma and urine. PCI found in blood originates from the liver and is capable of inhibiting several serine proteases involved in the regulation of coagulation and fibrinolysis,including activated protein C,thrombin,factor Xa,various kallikreins and plasminogen activators. PCI has been found to have antimicrobial and antitumor properties and thus appears to be a multi-functional protein.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF