



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human P-Selectin Monoclonal Antibody

单克隆抗体

产品货号: MA02100

第三版

免疫原:Fusion protein of P-selectin

亚型:IgG2b

克隆:756

描述:P-selectin,also known as GMP140 or CD62P,is a transmembrane glycoprotein that mediates the interaction of activated endothelial cells or platelets with leukocytes. It is an adhesion molecule involved in the pathogenesis of inflammation,thrombosis,and oncogenesis. P-selectin is stored in the alpha-granules of platelets and Weibel-Palade bodies of endothelial cells. Upon cell activation by agonists,P-selectin is transported rapidly to the cell surface.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF