



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human PSMD4 Monoclonal Antibody

单克隆抗体

产品货号: MA02101

第三版

免疫原:Fusion protein of PSMD4

亚型:IgG1

克隆:92

描述:SMD4,also named as MCB1,AF or ASF,belongs to the proteasome subunit S5A family. It binds and presumably selects ubiquitin-conjugates for destruction. PSMD4 displays selectivity for longer polyubiquitin chains. It modulates intestinal fluid secretion. PSMD4 is exposed on the sperm acrosomal surface,a notion that was further supported by the binding of non-cell permeant,biotinylated proteasomal inhibitor ZL3VS to the sperm acrosome. It is involved in the sperm-ZP penetration during fertilization. This antibody is specific to PSMD4.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF