



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RAB18 Monoclonal Antibody

单克隆抗体

产品货号: MA02103

第三版

免疫原:Fusion protein of Rab18

亚型:IgG1

克隆:89

描述:RAB18 is a member of a family of Ras-related small GTPases. RAB18 plays a role in apical endocytosis/recycling, and may be implicated in transport between the plasma membrane and early endosomes. It also plays a key role in eye and brain development and neurodegeneration. Defects in RAB18 are associated with Warburg micro syndrome type 3.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF