



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RAB27A Monoclonal Antibody

单克隆抗体

产品货号: MA02104

第三版

免疫原:Fusion protein of RAB27A

亚型:IgG1

克隆:200

**描述:**Rab27a is a member of a large family of Ras-related small GTPases. The protein is membrane-bound and may be involved in protein transport and small GTPase mediated signal transduction. Recent studies suggest that Rab27a is associated with melanosomes in pigmented cells and regulates melanosome transport via its interaction with actin-based cellular motors such as myosin Va and myosin VIIa. RAB27A is also required for regulated secretion in cytotoxic T lymphocytes. Rab27a is implicated in several diseases,including Griscelli syndrome,choroideremia,and the Hermansky-Pudlak syndrome mouse model,gunmetal.

**配方:**

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF