

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse RAC1 Monoclonal Antibody

单克隆抗体

产品货号: MA02106

第三版

免疫原:Fusion protein of Rac1

亚型:IgG2b

克隆:116

描述:RAC1b,also named as TC25,belongs to the small GTPase superfamily and Rho family. RAC1 is a plasma membrane-associated small GTPase which cycles between active GTP-bound and inactive GDP-bound states. In its active state,it binds to a variety of effector proteins to regulate cellular responses such as secretory processes,phagocytosis of apoptotic cells,epithelial cell polarization and growth-factor induced formation of membrane ruffles.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF