

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RBM4 Monoclonal Antibody

单克隆抗体

产品货号: MA02109

第三版

免疫原:Fusion protein of RBM4

亚型:IgG2a

克隆:634

描述:RBM4,RNA-binding motif protein4,contains 2 RRM-type RNA-binding motif and a retroviral-type(RT) zinc finger. RNA-binding factor participates in number of aspects of cellular processes such as alternative splicing of pre-mRNA and translation regulation. RBM4 recruits eIF4A1 to stimulate IRES-dependent translation in responds to cellular stress. Once assembled at the rHRE,the HIF2A-RBM4-eIF4E2 complex captures the 5-prime cap and targets mRNAs to polysomes for active translation,thereby evading hypoxia-induced repression of protein synthesis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF