



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse SCPEP1 Monoclonal Antibody

单克隆抗体

产品货号: MA02112

第三版

免疫原:Fusion protein of RISC

亚型:IgG2b

克隆:655

描述:Scpep1(Retinoid-inducible serine carboxypeptidase) is a 51 kDa protein which is also named as RISC and belongs to the peptidase S10 family. This protein may be involved in vascular wall and kidney homeostasis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF