



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human RNH1 Monoclonal Antibody

单克隆抗体

产品货号: MA02113

第三版

免疫原:Fusion protein of RNH1

亚型:IgG2b

克隆:347

描述:Ribonuclease/ angiogenin inhibitor (RNH,synonyms: RAI,MGC4569,MGC18200,MGC54054) is a protein that binds tightly to ribonucleases in cells and may be essential in the control of mRNA degradation and gene expression. The human RNH gene has been regionally localized to chromosome band 11p15 by in situ hybridization.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF