



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human S100A10 Monoclonal Antibody

单克隆抗体

产品货号: MA02117

第三版

免疫原:Fusion protein of S100A10

亚型:IgG1

克隆:379

描述:S100A10,also known as p11,is a member of the S100 family of small,EF hand containing dimeric proteins. S100 proteins are localized in the cytoplasm and/or nucleus of a wide range of cells,and involved in the regulation of a number of cellular processes such as cell cycle progression and differentiation. S100A10 is present on the surface of endothelial and other cells in a heterotetrameric complex with another Ca(2+)-binding protein,annexin II. S100A10 may function in exocytosis and endocytosis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF