



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse S100A11 Monoclonal Antibody

单克隆抗体

产品货号: MA02118

第三版

免疫原:Fusion protein of S100A11

亚型:IgG1

克隆:711

描述:S100A11,also named as MLN 70,Calgizzarin and S100C,Belongs to the S-100 family. It facilitates the differentiation and the cornification of keratinocytes. The naming of S100A11 systematically arises from its membership of the S100 protein family,named after their solubility in 100% saturated ammonium sulfate solution. First discovered in 1989,S100A11 has since been proposed to have direct biological functions in an assortment of physiological processes such as endo- and exocytosis,regulation of enzyme activity,cell growth and regulation,apoptosis and inflammation.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF