

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human S100A6 Monoclonal Antibody

单克隆抗体

产品货号: MA02119

第三版

免疫原:Fusion protein of S100A6

亚型:IgG1

克隆:127

描述:S100A6 is also named as calyculin, prolactin receptor-associated protein (PRA), growth factor-inducible protein 2A9 or MLN4. It belongs to S100 family of low molecular weight, acidic, calcium-binding proteins which contain two EF-hand calcium binding sites. S100A6 may function as a calcium sensor to activate several processes in the calcium signal transduction pathway of cell growth, proliferation, secretion and exocytosis.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF