



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SECTM1 Monoclonal Antibody

单克隆抗体

产品货号: MA02122

第三版

免疫原:Fusion protein of SECTM1

亚型:IgM

克隆:501

描述:SECTM1,also named as protein K12,is originally identified by its location just upstream of the CD7 locus. It has been proposed as a ligand for CD7. SECTM1 is involved in thymocyte signaling. This antibody is specific to SECTM1.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF