

---

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

---

Anti-Human/Mouse SEPT2 Monoclonal Antibody

单克隆抗体

产品货号: MA02123

---

第三版

免疫原:Fusion protein of SEPT2

亚型:IgG1

克隆:260

描述:SEPT2,also named as DIFF6,KIAA0158,NEDD5,belongs to the septin family. It is a filament-forming cytoskeletal GTPase. SEPT2 is required for normal organization of the actin cytoskeleton. It plays a role in ciliogenesis and collective cell movements. SEPT2 is coordinated expression with SEPT6 and SEPT7.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF