



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SEPT8 Monoclonal Antibody

单克隆抗体

产品货号: MA02124

第三版

免疫原:Fusion protein of SEPT8

亚型:IgM

克隆:1155

描述:SEPT8,also named as KIAA0202,is a member of the highly conserved septin family. Septins are membrane-associated GTPases which are involved in cytokinesis and cellular morphogenesis. SEPT8 is an interaction partner of SEPT4. It may play a role in platelet secretion.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF