

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human APCS Monoclonal Antibody

单克隆抗体

产品货号: MA02126

第三版

免疫原:Fusion protein of Serum amyloid P component

亚型:IgG2b

克隆:4655

描述:Serum amyloid P component (SAP),a member of pentraxin family,is a constitutive serum protein that is synthesized by hepatocytes. It is an acute phase protein,structurally related to C-reactive protein. SAP is a calcium-dependent ligand binding protein,which can interact with DNA and histones and may scavenge nuclear material released from damaged circulating cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF