

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SLAMF1 Monoclonal Antibody

单克隆抗体

产品货号: MA02129

第三版

免疫原:Fusion protein of SLAM/CD150

亚型:IgM

克隆:834

描述:Signaling lymphocytic activation molecule (SLAM,also known as SLAMF1 or CD150) is a 70-kDa transmembrane glycoprotein that belongs to the immunoglobulin gene superfamily. SLAM is constitutively expressed on peripheral blood memory T-cells,T-cell clones,immature thymocytes and a proportion of B-cells,and is rapidly induced on naive T-cells after activation. SLAM is involved in T cell stimulation and cytokine production,and may play an important role in the regulation of immune responses to pathogens.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF