

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SMARCD1 Monoclonal Antibody

单克隆抗体

产品货号: MA02130

第三版

免疫原:Fusion protein of SMARCD1

亚型:IgG1

克隆:224

描述:SMARCD1 (SWI/SNF related,matrix associated,actin dependent regulator of chromatin,subfamily d,member 1),also known as BAF60a (BRG1-associated factor 60A),Rsc6p or CRACD1,is a 476 amino acid protein that localizes to the nucleus and contains one SWIB domain. Expressed ubiquitously with notable expression in liver,brain,muscle, lung,kidney,pancreas and placenta,SMARCD1 functions as a component of the BAF complex and,in conjunction with a variety of other proteins,plays an essential role in chromatin remodeling. In addition,SMARCD1 influences vitamin D-mediated transcriptional activity and is thought to provide a link between the vitamin D receptor (VDR) and SWI/SNF chromatin remodeling complexes. Multiple isoforms of BAF60a exist due to alternative splicing events.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF