

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat SMN1 Monoclonal Antibody

单克隆抗体

产品货号: MA02132

第三版

免疫原:Fusion protein of SMN (Human,Mouse,Rat)

亚型:IgG2b

克隆:600

描述:The survival of motor neurons (SMN) genes are the disease genes of spinal muscular atrophy (SMA),a common motor neuron degenerative disease. The level of SMN protein correlates with phenotypic severity of SMA. SMA patients lack a functional SMN1 gene,but they possess an intact SMN2 gene,which though nearly identical to SMN1,is only partially functional,because a large majority of SMN2 transcripts lack exon 7,resulting in production of a truncated,less stable SMN protein. This antibody can recognize human,mouse and rat SMN gene.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF