



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SNRPD2 Monoclonal Antibody

单克隆抗体

产品货号: MA02137

第三版

免疫原:Fusion protein of SNRPD2

亚型:IgG2b

克隆:9321

描述: Systemic lupus erythematosus is characterized by antibodies to a variety of intracellular self-antigens, such as dsDNA and Sm, and these serve as hallmarks in the diagnosis of systemic autoimmune diseases. SNRPD2 is one of Sm protein, and is required for pre-mRNA splicing, snRNP biogenesis.

配方:

浓度: None

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF