

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human SRP9 Monoclonal Antibody

单克隆抗体

产品货号: MA02138

第三版

免疫原:Fusion protein of SRP9

亚型:IgG2a

克隆:124

描述:Signal recognition particle 9(SRP9) is component of the signal recognition particle (SRP),which is a ribonucleoprotein complex that mediates the targeting of proteins to the rough endoplasmic reticulum (ER). SRP9 form a heterodimer with SRP14,which involves in arrest of the elongation of the nescent chains during targeting to ensure efficient translocation of the preprotein.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF