

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-Human/Mouse/Rat ST2 Monoclonal Antibody

单克隆抗体

产品货号: MA02140

第三版

免疫原:Fusion protein of ST2

亚型:IgG1

克隆:6478

描述:Protein ST2,also known as IL1RL1, is a member of the interleukin 1 receptor family. The gene of ST2 encodes three splice variants: a soluble secreted form,a transmembrane form and a form of undetermined localization.This antibody recognizes at least 80kd,40kd and 30kd forms in different tissues and cells.

配方:

浓度:None

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF